

Catheter Carts



Specifications

Bulk Storage Catheter Carts

ITEM # _____

PROJECT _____

For catheter storage in the E.R., O.R. or Central Supply, this versatile cart can hold catheters up to 60" long. This cart comes in 3 sizes and allows for both hanging storage as well as shelf storage. Label holders, bins and ledges come standard.


- Overall height: 69"
- 5" swivel caster with polyurethane tread, 2 with total-lock brakes
- Zinc or Chrome finish
- Cart covers available in uncoated, coated and antimicrobial material
- Ships assembled or K.D.



CPC-36-BS



CPC-48-BS



CPC-60-BS

Catheter Carts

W x L (in.)	W x L (mm)	Style	# of Bars per Rack	# of Hooks Included	Max. Hooks per Bar	# of Totes	# of Bins	# of Ledges	# of Label Holders	Lacquered Zinc Model #	Chrome Model #
24" x 36"	610 x 1524	Bulk Storage	4 – 36" Bars	24	9	-	4	6	4	CPC-36LZ-BS	CPC-36-BS
24" x 48"	610 x 1219	Bulk Storage	4 – 48" Bars	32	13	-	6	6	6	CPC-48LZ-BS	CPC-48-BS
24" x 60"	610 x 1524	Bulk Storage	4 – 60" Bars	40	17	-	6	6	6	CPC-60LZ-BS	CPC-60-BS



IRSG
12230 NE Woodinville Drive, Suite A
Woodinville, WA 98072
Phone: (425) 398-1875 • Fax: (425) 398-1675

IRSG

For product information, call 1-888-398-2009

For custom design, configuration or fabrication, please call (888) 398-2009 or e-mail us at info.irsg.com
For more information please visit www.irsg.com

Every attempt has been made to ensure the accuracy of the product information. We cannot be held responsible for printing or typographical errors. Please confirm product information a time of order; information is subject to change without notice