

Catheter Carts



Specifications

*Combination
Catheter Carts - Economy*

ITEM # _____

PROJECT _____



**Economy Catheter Cart
CPCC-60-24E
Shown with optional
Catheter Guard Rail**



**Economy Catheter Cart
CPCC-60-36E
Shown with optional
Catheter Guard Rail**

Catheter Carts

W x L (in.)	W x L (mm)	Style	# of Bars per Rack	# of Hooks Included	Max. Hooks per Bar	# of Totes	# of Slides	# of Bins	# of Ledges	# of Label Holders	Catheter Guard Rail*	Lacquered Zinc Model #	Chrome Model #
24" x 60"	610 x 1524	3 Side Shelves, 4 - 24" Bars	4	16	5	-	-	-	6	-	CGR-24D	CPCC-60LZ-24E	CPCC-60-24E
24" x 60"	610 x 1524	3 Side Shelves, 4 - 36" Bars	4	24	9	-	-	-	6	-	CGR-36D	CPCC-60LZ-36E	CPCC-60-36E

* Overall Height - 69"



IRSG

IRSG
12230 NE Woodinville Drive, Suite A
Woodinville, WA 98072
Phone: (425) 398-1875 • Fax: (425) 398-1675

For product information, call 1-888-398-2009

For custom design, configuration or
fabrication, please call (888) 398-2009 or
e-mail us at info.irsg.com
For more information
please visit www.irsg.com

Every attempt has been made to ensure the accuracy of the product information. We cannot be held responsible for printing or typographical errors. Please confirm product information a time of order; information is subject to change without notice