

Catheter Carts



Specifications

Combination Catheter Carts

- Standard or heavy-duty models available
- Heavy-duty models come standard with 6" casters, 2 with brakes
- Standard models come standard with 5" polyurethane swivel casters, 2 with total-lock brakes
- 3" non-marring bumpers

- All components are adjustable
- Clear plastic label holders attach to end of hook to identify catheters
- Optional catheter guard *
- Ships assembled or K.D. (easy assembly)

ITEM # _____

PROJECT _____



CPC-60-36
Standard Combination
Catheter Cart



CPC-60-D36
Heavy-Duty
Deluxe Catheter Cart

Catheter Carts

W x L (in.)	W x L (mm)	Style	# of Bars per Rack	# of Hooks Included	Max. Hooks per Bar	# of Totes	# of Tote Slides	# of Bins	# of Ledges	# of Label Holders	Catheter Guard Rail*	Lacquered Zinc Model #	Chrome Model #
24" x 48"	610 x 1219	3 Side Shelves, 4 – 24" Bars	4	16	5	2	1	6	6	6	CGR-48	CPCC-48-24LZ	CPCC-48-24
24" x 60"	610 x 1524	3 Side Shelves, 4 – 24" Bars	4	16	5	2	1	6	6	6	CGR-60	CPCC-60-24LZ	CPCC-60-24
24" x 60"	610 x 1524	3 Side Shelves, 4 – 36" Bars	4	24	9	2	1	6	6	6	CGR-60	CPCC-60-36LZ	CPCC-60-36

Heavy - Duty Model:

W x L (in.)	W x L (mm)	Style	# of Bars per Rack	# of Hooks Included	Max. Hooks per Bar	# of Totes	# of Tote Slides	# of Bins	# of Ledges	# of Label Holders	Catheter Guard Rail*	Lacquered Zinc Model #	Chrome Model #
24" x 60"	610 x 1524	3 Side Shelves, 4 – 24" Bars	4	16	5	2	1	6	6	6	CGR-24D	CPCC-60-D24LZ	CPCC-60-D24
24" x 60"	610 x 1524	3 Side Shelves, 4 – 36" Bars	4	24	9	2	1	6	6	6	CGR-36D	CPCC-60-D36LZ	CPCC-60-D36



IRSG
12230 NE Woodinville Drive, Suite A
Woodinville, WA 98072
Phone: (425) 398-1875 • Fax: (425) 398-1675

IRSG

For product information, call 1-888-398-2009

For custom design, configuration or fabrication, please call (888) 398-2009 or e-mail us at info.irsg.com
For more information please visit www.irsg.com

Every attempt has been made to ensure the accuracy of the product information. We cannot be held responsible for printing or typographical errors. Please confirm product information a time of order; information is subject to change without notice