

Inhalation Therapy Cart



Specifications

Inhalation Therapy Cart

ITEM # _____

PROJECT _____



- Construction – Quick assembly without the use of tools. Shelf heights are adjustable in 1" (25.4mm) increments.
- Finish – Handles, frame and bottom shelf are chrome-plated. Top and bottom inserts are PVC-coated to reduce noise and vibration.
- Size - Divider grid can accommodate twelve D & E size cylinders. Each PVC-coated compartment measures 4 1/2" (114mm) square.
- Casters - Chrome-plated friction stem style - 5" (127mm) conductive casters, two total lock and two free swivels.
- Design – Conductive split-sleeves supports help reduce static electricity build up.
- Ships K.D.

W x L (in.)	W x L (mm)	Height (in.) (mm)	Lacquered Zinc Model #	Chrome Model #
18" x 24"	457 x 610	40" 1016	GCC-1824LZ	GCC-1824C

Inhalation Therapy Cart



IRSG
12230 NE Woodinville Drive, Suite A
Woodinville, WA 98072
Phone: (425) 398-1875 • Fax: (425) 398-1675

IRSG For product information, call 1-888-398-2009

For custom design, configuration or fabrication, please call (888) 398-2009 or e-mail us at info.irsg.com
For more information please visit www.irsg.com