

Oxygen Cart



Specifications *Oxygen Cart*

ITEM # _____

PROJECT _____

Portable, lightweight compressed oxygen system featuring a rubber handgrip for handling comfort and a locking bracket for cylinder stability.



- Lightweight steel tubing cart with chrome finish
- Ergonomic height adjustments from 37" to 24 ½"
- Accommodates D or E cylinders
- 2 - 6" non-marring tires
- Support stabilizer with two reinforced rubber tips prevents tipping over when stationary
- Perfect for both indoor and outdoor use
- Available in Zinc or Chrome

W x L x H (in.)	W x L x H (mm)	Model #
4 x 22 x 4	102 x 559 x 102	OTC

OTC - Cylinder not included

Oxygen Cart



IRSG
12230 NE Woodinville Drive, Suite A
Woodinville, WA 98072
Phone: (425) 398-1875 • Fax: (425) 398-1675

IRSG For product information, call 1-888-398-2009

For custom design, configuration or fabrication, please call (888) 398-2009 or e-mail us at info.irsg.com
For more information please visit www.irsg.com

Every attempt has been made to ensure the accuracy of the product information. We cannot be held responsible for printing or typographical errors. Please confirm product information a time of order; information is subject to change without notice