

Patient Care Cart



Specifications

Patient Care Cart



PCCH-2448C Heavy-Duty Patient Care Cart



PCC-2448C Standard Patient Care Cart

ITEM # _____

PROJECT _____

- 6 shelves, 5 featuring 11 - 6" high adjustable shelf dividers and an open bottom shelf
- Shelves, posts and dividers have a chrome-plated finish
- Open wire construction reduces dirt and dust build-up.
- 5" non-marring casters, two swivel and two with total-lock brakes and polyurethane tread
- Quick assembly - shelf heights are adjustable in 1" increments
- 1,000 lb. total load capacity

Patient Care Cart

W x L (In.)	W x L (mm)	Overall Height (In.)	Overall Height (mm)	Style	Hang Rod Length	Shelf Length	Label Holders	Enclosure Rods	Dividers	Lacquered Zinc Model #	Chrome Model #
24" x 48"	610 x 1219	69"	1752	Standard	24	24	8	10	8	PCC-2448LZ	PCC-2448C
24" x 60"	610 x 1524	69"	1752	Standard	24	36	12	12	12	PCC-2460LZ	PCC-2460C
24" x 48"	610 x 1219	71.5"	1861	Heavy-Duty	24	24	8	10	8	PCCH-2448LZ*	PCCH-2448C*
24" x 60"	610 x 1524	71.5"	1861	Heavy-Duty	24	36	12	12	12	PCCH-2460LZ*	PCCH-2460C*

*Dollies add 4" to length and 2" to width to accommodate the continuous bumpers.



IRSG

IRSG
12230 NE Woodinville Drive, Suite A
Woodinville, WA 98072
Phone: (425) 398-1875 • Fax: (425) 398-1675

For product information, call 1-888-398-2009

For custom design, configuration or fabrication, please call (888) 398-2009 or e-mail us at info.irsg.com
For more information please visit www.irsg.com

Every attempt has been made to ensure the accuracy of the product information. We cannot be held responsible for printing or typographical errors. Please confirm product information a time of order; information is subject to change without notice