

# Peel Pack Cart



## Specifications *Peel Pack Cart*

ITEM # \_\_\_\_\_

PROJECT \_\_\_\_\_



- 6 shelves, 5 featuring 11 - 6" high adjustable shelf dividers and an open bottom shelf
- Shelves, posts and dividers have a chrome-plated finish
- Open wire construction reduces dirt and dust build-up.
- 5" non-marring casters, two swivel and two with total-lock brakes and polyurethane tread
- Quick assembly - shelf heights are adjustable in 1" increments
- 1,000 lb. total load capacity

Peel Pack Cart

Width x Length	Overall Height	Model No.
24 x 60	69"	PPC-2460C



IRSG  
12230 NE Woodinville Drive, Suite A  
Woodinville, WA 98072  
Phone: (425) 398-1875 • Fax: (425) 398-1675

**IRSG** For product information, call 1-888-398-2009

For custom design, configuration or fabrication, please call (888) 398-2009 or e-mail us at [info.irsg.com](mailto:info.irsg.com)  
For more information please visit [www.irsg.com](http://www.irsg.com)